

# Table Of Contents

<b>C<sub>18</sub></b> .....	<b>Pages 3-43</b>
500Å .....	Pages 3-43
<b>DVB</b> .....	<b>Pages 45-105</b>
500Å .....	Pages 45-91
10 <sup>3</sup> Å.....	Pages 92-103
Solid Bead.....	Pages 104-105
<b>Glucose DVB</b> .....	<b>Pages 107-108</b>
500Å.....	Pages 107-108
<b>Hydroxylate DVB</b> .....	<b>Pages 110</b>
500Å.....	Pages 110
<b>Organic Acid DVB</b> .....	<b>Pages 112</b>
500Å.....	Pages 112
<b>Polyamino DVB</b> .....	<b>Pages 114-118</b>
500Å.....	Pages 114-118
<b>Peptide Protein DVB</b> .....	<b>Pages 120-132</b>
10 <sup>3</sup> Å.....	Pages 120
10 <sup>4</sup> Å.....	Page 121-132
<b>Sulfonated DVB</b> .....	<b>Pages 134-136</b>
500Å .....	Page 134-136
<b>Sax Quat DVB</b> .....	<b>Pages 138-141</b>
10 <sup>3</sup> Å.....	Page 138-141
<b>Index</b> .....	<b>Pages 142-145</b>



**C<sub>18</sub> DVB**






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## WATER SOLUBLE VITAMINS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 50/40/10 H<sub>2</sub>O/ACN/MeOH  
 pH 4.0 w/ CH<sub>3</sub>COOH  
**Flow Rate:** 0.7mL/min.  
**Injection:** 15µL  
**Temperature:** 25°C  
**Detector:** UV @246nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## WATER SOLUBLE VITAMINS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 50/40/10 H<sub>2</sub>O/ACN/MeOH  
pH 4.0 w/ CH<sub>3</sub>COOH  
**Flow Rate:** 0.7mL/min.  
**Injection:** 15µL  
**Temperature:** 25°C  
**Detector:** UV @206nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## WATER SOLUBLE VITAMINS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 70/30 H<sub>2</sub>O/ACN w/ 0.1% TFA  
**Flow Rate:** 0.7mL/min.  
**Injection:** 20µL  
**Temperature:** 25°C  
**Detector:** UV @268nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## WATER SOLUBLE VITAMINS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 50/40/10 H<sub>2</sub>O/ACN/MeOH  
pH 4.0 w/CH<sub>3</sub>COOH  
**Flow Rate:** 0.7mL/min.  
**Injection:** 20µL  
**Temperature:** 25°C  
**Detector:** UV @245nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## FAT SOLUBLE VITAMINS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** ACN w/ 0.1% TFA  
**Flow Rate:** 1.5mL/min.  
**Injection:** 50µL  
**Temperature:** 25°C  
**Detector:** UV @284nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

ASPIRIN and RELATED COMPOUNDS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 60/40 ACN/H<sub>2</sub>O w/ 0.1% TFA  
**Flow Rate:** 0.7mL/min.  
**Injection:** 2µL  
**Temperature:** 25°C  
**Detector:** UV @235nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## BUTYL PHENOL STANDARDS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 70/30 ACN/H<sub>2</sub>O w/ 0.1% TFA  
**Flow Rate:** 0.7mL/min.  
**Injection:** 5µL  
**Temperature:** 25°C  
**Detector:** UV @215nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## CHOLINE COMPOUNDS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 90/3/7 0.1M Na<sub>2</sub>HPO<sub>4</sub>/ACN/MeOH  
**Flow Rate:** 0.7mL/min.  
**Injection:** 10µL  
**Temperature:** 25°C  
**Detector:** UV @215nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## COLUMN TEST MIX

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 65/30/5 ACN/H<sub>2</sub>O/MeOH w/0.1%TFA  
**Flow Rate:** 0.7mL/min.  
**Injection:** 5µL  
**Temperature:** 25°C  
**Detector:** UV @254nm







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## VASOCONSTRICTORS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 24/75/1 ACN/0.2M NaOH/Butylamine  
**Flow Rate:** 0.7mL/min.  
**Injection:** 10µL  
**Temperature:** 25°C  
**Detector:** UV @254nm







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## VANILLIN COMPOUNDS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 40/60 H<sub>2</sub>O/ACN w/0.1% TFA  
**Flow Rate:** 0.7mL/min.  
**Injection:** 4µL  
**Temperature:** 25°C  
**Detector:** UV @235nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## NUCLEOSIDES and BASES

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 99/1 0.01%M Sodium Acetate/ACN  
pH 4.0 w/TFA  
**Flow Rate:** 0.7mL/min.  
**Injection:** 15µL  
**Temperature:** 25°C  
**Detector:** UV @254nm







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

NUCLEOSIDES and BASES  
(at different pH's)

Part Number: 18500  
 Packing: Jordi DVB C<sub>18</sub> 500Å  
 Column: 15cm X 4.6mm ID  
 Solvent: 99/1 0.01%M Sodium Acetate/ACN  
 pH 4.0 w/TFA & pH 7.1 w/TFA  
 Flow Rate: 0.7mL/min.  
 Injection: 15µL  
 Temperature: 25°C  
 Detector: UV @254nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

NUCLEOSIDES and BASES  
(at different Temperatures)

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 99/1 0.01%M Sodium Acetate/ACN  
 pH 4.0 w/TFA  
**Flow Rate:** 0.7mL/min.  
**Injection:** 15µL  
**Temperature:** See curve detail  
**Detector:** UV @254nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ANTICONVULSANTS

**Part Number:** 18500  
**Packing:** DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent A:** H<sub>2</sub>O  
**Solvent B:** 87.5/12.5 ACN/MeOH  
**Flow Rate:** 0.7mL/min.  
**Injection:** 20µL  
**Temperature:** 25°C  
**Detector:** UV @215nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## BARBITURATES

**Part Number:** 18500  
**Packing:** DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent :** 70/22.5/7.5 0.05 M Na<sub>2</sub>HPO<sub>4</sub>/ACN/MeOH  
**Flow Rate:** 0.7mL/min.  
**Injection:** 12µL  
**Temperature:** 25°C  
**Detector:** UV @210nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## CHOLINE COMPOUNDS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** H<sub>2</sub>O @pH 3.3 w/HCl  
**Flow Rate:** 0.7mL/min.  
**Injection:** 15µL  
**Temperature:** 25°C  
**Detector:** UV @200nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## VANILLIN COMPOUNDS

**Part Number:** 18501  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 40/60 H<sub>2</sub>O/ACN w/0.1% TFA  
**Flow Rate:** 0.7mL/min.  
**Injection:** 4µL  
**Temperature:** 25°C  
**Detector:** UV @235nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ASPIRIN and RELATED COMPOUNDS

**Part Number:** 18501  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 40/60 H<sub>2</sub>O/ACN w/0.1% TFA  
**Flow Rate:** 0.7mL/min.  
**Injection:** 2µL  
**Temperature:** 25°C  
**Detector:** UV @235nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.


# HPLC APPLICATION

## 2' and 3' MONOPHOSPHATE NUCLEOTIDES

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 99/1 0.01M Sodium Acetate/ACN  
 w/0.1% TFA  
**Flow Rate:** 0.7mL/min.  
**Injection:** 10µL  
**Temperature:** 25°C  
**Detector:** UV @254nm



- 1. 3'-AMP
- 2. 2'-AMP
- 3. 2'-CMP & 3'-CMP
- 4. 3'-AMP
- 5. 3'-UMP
- 6. 2'-AMP & 2' -UMP
- 7. 3' -GMP
- 8. 2' -GMP







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## 2' and 3' MONOPHOSPHATE NUCLEOTIDES

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 99/1 0.01M Sodium Acetate/ACN  
w/0.1% TFA  
**Flow Rate:** 0.7mL/min.  
**Injection:** 7µL  
**Temperature:** 25°C  
**Detector:** UV @254nm



1. 3'-UMP
2. 2'-UMP
3. 3'-GMP
4. 2'-GMP






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PHENOLS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent A:** 44/56 H<sub>2</sub>O/ACN w/0.1% TFA  
**Solvent B:** 10/90 H<sub>2</sub>O/ACN w/0.1% TFA  
**Flow Rate:** 0.7mL/min.  
**Injection:** 2µL of 2mg/mL solution  
**Temperature:** 25°C  
**Detector:** UV @280nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## CATECHOLAMINES

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent A:** 98/1/1 0.01M NaH<sub>2</sub>PO<sub>4</sub>/ACN/Butylamine  
**Solvent B:** 79/20/1 0.01M NaH<sub>2</sub>PO<sub>4</sub>/ACN/Butylamine  
**Flow Rate:** 0.7mL/min.  
**Injection:** 50µL  
**Temperature:** 25°C  
**Detector:** UV @280nm







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ALKALOID DRUGS

**Part Number:** 18503  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 10cm X 4.6mm ID  
**Solvent:** 70/20/10 H<sub>2</sub>O/ACN/MeOH w/0.1% TFA  
**Flow Rate:** 0.5mL/min.  
**Injection:** 25µL  
**Temperature:** 25°C  
**Detector:** UV @230nm



1. Codeine (0.14 mg/mL)
2. Oxycodone (0.29 mg/mL)
3. Meperidine (0.57 mg/mL)




MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## STERIODS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 65/35 ACN/H<sub>2</sub>O  
**Flow Rate:** 0.7mL/min.  
**Injection:** 20µL  
**Temperature:** 25°C  
**Detector:** UV @215nm, 0.62 AUFS






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## SHORT CHAIN FATTY ACIDS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 50/50 ACN/H<sub>2</sub>O/H<sub>3</sub>PO<sub>4</sub>  
**Flow Rate:** 0.7mL/min.  
**Injection:** 150µL  
**Temperature:** 25°C  
**Detector:** UV @210nm, 0.6 AUFS








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## TETRACYCLINE and RELATED COMPOUNDS

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 20/80 ACN/H<sub>2</sub>O w/0.1% TFA  
**Flow Rate:** 0.7mL/min.  
**Injection:** 20µL  
**Temperature:** 25°C  
**Detector:** UV @300nm, 0.5 AUFS







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ERYTHROMYCIN (New and Aged)

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 45/23/32 ACN/0.05M KH<sub>2</sub>PO<sub>4</sub>/MeOH  
**Flow Rate:** 0.7mL/min.  
**Injection:** 35µL of 6mg/mL  
**Temperature:** 25°C  
**Detector:** RI (Waters 401) @8X







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

PENICILLIN  
(New and Aged)

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 50/50 ACN/0.05M KH<sub>2</sub>PO<sub>4</sub>  
**Flow Rate:** 0.7mL/min.  
**Injection:** 1.5µL of 17mg/mL  
**Temperature:** 25°C  
**Detector:** UV @220nm 2.0 AUFS








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## AMPICILLIN

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 50/50 ACN/0.05M KH<sub>2</sub>PO<sub>4</sub>  
**Flow Rate:** 0.7mL/min.  
**Injection:** 3μL of 14mg/mL  
**Temperature:** 25°C  
**Detector:** UV @210nm 2.0 AUFS







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

FUMONISIN B<sub>1</sub>  
(From *Fusarium Moniliforme*)

**Part Number:** 18501  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 50/25/25 H<sub>2</sub>O/ACN/MeOH  
 pH=3.5 w/ Acetic Acid  
**Flow Rate:** 1.0mL/min.  
**Injection:** 250µL of 250ppm Solution in  
 MeOH/H<sub>2</sub>O (80/20)  
**Temperature:** 40°C  
**Detector:** Evaporative Light Scattering







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

**FUMONISIN B<sub>1</sub>**  
(From *Fusarium Moniliforme*)

**Part Number:** 18501  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 50/25/25 H<sub>2</sub>O/MeOH/ACN  
 pH=3.0 w/Acetic Acid  
**Flow Rate:** 1.0mL/min.  
**Injection:** 250µL of 250ppm Solution in  
 MeOH/H<sub>2</sub>O (80/20)  
**Temperature:** 40°C  
**Detector:** Evaporative Light Scattering



ELSD MK III  
 Tube T=100 deg C  
 Gas=2.0 SLPM  
 ATTN=16  
 Time Const.=1 Sec







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

**FUMONISIN B<sub>1</sub>**  
(From *Fusarium Moniliforme*)

**Part Number:** 18501  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 50/25/25 H<sub>2</sub>O/MeOH/ACN  
 pH=4.0/5.0 w/Acetic Acid  
**Flow Rate:** 1.0mL/min.  
**Injection:** 250µL of 250ppm Solution in MeOH/  
 H<sub>2</sub>O (80/20)  
**Temperature:** 40°C  
**Detector:** Evaporative Light Scattering



ELSD MK III  
 Tube T=100 deg C  
 Gas=2.0 SLPM  
 ATTN=16  
 Time Const.=1 Sec







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## REVERSE PHASE STANDARD MIX (New Mobile Phase)

**Part Number:** 18501  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 60/40 THF/H<sub>2</sub>O  
**Flow Rate:** 1.0mL/min.  
**Injection:** 5µL  
**Temperature:** 50°C  
**Detector:** UV @254nm



Note: The importance of solvent/column interaction using Jordi DVB columns *cannot* be overemphasized. We have found that a 50/50/ mix of MeOH/ACN for the strong solvent is adequate for many Reverse-Phase separations and is better than either alone. We have now observed that the use of THF/ACN as strong solvent is often better than MeOH/ACN. In general Lewis bases (electron donor solvents) deactivate the aromatic rings and often dramatically increase the column efficiencies. The following data describes some Reverse Phase results for our 500Å DVB Column.




MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ASCORBIC and DEHYDROASCORBIC ACID

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 40:60 Water (pH 2.7 w/ Acetic Acid)/Acetonitrile/Methanol (5:1)  
**Flow Rate:** 0.7mL/min.  
**Injection:** 10µL  
**Concentration:** 1mg/mL  
**Temperature:** Ambient  
**Detector:** ELSD









MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## RP MIX

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 10/20/30/40 MeOH/ACN/THF/H<sub>2</sub>O  
**Flow Rate:** 1.0mL/min.  
**Injection:** 5µL  
**Concentration:** N/A  
**Temperature:** 50°C  
**Detector:** UV @254nm, Sens. 2.0 AUFS, Press. 1900 psig






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.




# HPLC APPLICATION

## REVERSE PHASE STANDARD MIX (New Mobile Phase)

**Part Number:** 18501  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 50/25/25 ACN/MeOH/H<sub>2</sub>O  
**Flow Rate:** 1.0mL/min.  
**Injection:** 1µL  
**Temperature:** Ambient  
**Detector:** UV @254nm



Note: The importance of solvent/column interaction using Jordi DVB columns *cannot* be overemphasized. We have found that a 50/50/ mix of MeOH/ACN for the strong solvent is adequate for many Reverse-Phase separations and is better than either alone. We have now observed that the use of THF/ACN as strong solvent is often better than MeOH/ACN. In general Lewis bases (electron donor solvents) deactivate the aromatic rings and often dramatically increase the column efficiencies. The following data describes some Reverse Phase results for our 500Å DVB Column.









MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## REVERSE PHASE STANDARD MIX (New Mobile Phase)

**Part Number:** 18501  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 30/30/40 THF/ACN/H<sub>2</sub>O  
**Flow Rate:** 1.0mL/min.  
**Injection:** 5µL  
**Temperature:** 60°C  
**Detector:** UV @254nm



Note: The importance of solvent/column interaction using Jordi DVB columns *cannot* be overemphasized. We have found that a 50/50/ mix of MeOH/ACN for the strong solvent is adequate for many Reverse-Phase separations and is better than either alone. We have now observed that the use of THF/ACN as strong solvent is often better than MeOH/ACN. In general Lewis bases (electron donor solvents) deactivate the aromatic rings and often dramatically increase the column efficiencies. The following data describes some Reverse Phase results for our 500Å DVB Column.






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## REVERSE PHASE MIX

**Part Number:** 18500  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 15cm X 4.6mm ID  
**Mobile Phase:** 10/20/30/40 MeOH/ACN/THF/H<sub>2</sub>O  
**Flow Rate:** 1.0mL/min.  
**Injection:** 5µL  
**Temperature:** 50°C  
**Detector:** UV @254nm, 2.0 AUFS







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

FUMONISIN B<sub>1</sub>  
(From Fusarium Moniliforme)

**Part Number:** 18501  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 50/25/25 H<sub>2</sub>O/MeOH/ACN  
 pH=2.7 w/Formic Acid  
**Flow Rate:** 1.0mL/min.  
**Injection:** 250µL of 250ppm Solution in  
 MeOH/H<sub>2</sub>O (80/20)  
**Temperature:** 40°C  
**Detector:** Evaporative Light Scattering,






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## REVERSE PHASE STANDARD MIX (New Mobile Phase)

**Part Number:** 18501  
**Packing:** Jordi DVB C<sub>18</sub> 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 30/30/40 THF/ACN/H<sub>2</sub>O  
**Flow Rate:** 1.0mL/min.  
**Injection:** 5µL  
**Temperature:** 60°C  
**Detector:** UV @254nm





**DVB**










MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## CAFFEINE in COFFEE

**Part Number:** 16502  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 50/15/35 0.01M LiNO<sub>3</sub>/ACN/MeOH  
**Flow Rate:** 0.5mL/min.  
**Injection:** 50µL  
**Temperature:** 25°C  
**Detector:** UV @254nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## CAFFEINE in TEA & CHOCOLATE

**Part Number:** 16502  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 50/15/35 H<sub>2</sub>O/ACN/MeOH  
 H<sub>2</sub>O, 0.01M LiNO<sub>3</sub>  
**Flow Rate:** 0.5mL/min.  
**Injection:** 50µL  
**Temperature:** 50°C  
**Detector:** UV @254nm












MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## CAFFEINE in COKE & PEPSI

**Part Number:** 16502  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 15cm X 4.6mm ID  
**Solvent:** 50/15/35 0.01M LiNO<sub>3</sub>/ACN/MeOH  
**Flow Rate:** 0.5mL/min.  
**Injection:** 50µL  
**Temperature:** 50°C  
**Detector:** UV @254nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## BASE-NEUTRAL-ACID TEST MIX

**Part Number:** 16502  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 55/37.2/7.8 ACN/H<sub>2</sub>O/MeOH  
pH 3.0 w/0.05M Na<sub>2</sub>HPO<sub>4</sub>  
**Flow Rate:** 2.0mL/min.  
**Injection:** 40µL  
**Temperature:** 25°C  
**Detector:** UV @254nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## BUTYLPHENOL STANDARDS

**Part Number:** 16502  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 80/20 ACN/H<sub>2</sub>O w/0.05% TFA  
**Flow Rate:** 2.0mL/min.  
**Injection:** 20µL  
**Temperature:** 25°C  
**Detector:** UV @254nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ORGANIC ACIDS

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 90/10 MeOH/THF  
**Flow Rate:** 2.0mL/min.  
**Injection:** 50µL  
**Temperature:** 25°C  
**Detector:** UV @210nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ORGANIC ACIDS

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 34.7/15.3/50 ACN/THF/H<sub>2</sub>O  
pH 3.0 w/ H<sub>3</sub>PO<sub>4</sub>  
**Flow Rate:** 2.0mL/min.  
**Injection:** 150µL  
**Temperature:** 25°C  
**Detector:** UV @210nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## VASOCONSTRICTORS

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 24/74/2 ACN/0.2M NaOH/Butylamine  
**Flow Rate:** 3.0mL/min.  
**Injection:** 100µL  
**Temperature:** 25°C  
**Detector:** UV @254nm

1. Phenylephrine Hydrochloride
2. Norephedrine Hydrochloride
3. Procainamide Hydrochloride
4. Ephedrine Hydrochloride
5. Pseudoephedrine Hydrochloride








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ASPIRIN and RELATED COMPOUNDS

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 50/50 ACN/H<sub>2</sub>O w/0.05% TFA  
**Flow Rate:** 3.0mL/min.  
**Injection:** 100µL  
**Temperature:** 25°C  
**Detector:** UV @254nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## COLUMN TEST MIX

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 75/20/5 ACN/H<sub>2</sub>O/MeOH w/0.1%TFA  
**Flow Rate:** 2.0mL/min.  
**Injection:** 20µL  
**Temperature:** 25°C  
**Detector:** UV @254nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## VANILLIN COMPOUNDS

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 80/20 ACN/H<sub>2</sub>O w/0.1%TFA  
**Flow Rate:** 3.0mL/min.  
**Injection:** 10µL  
**Temperature:** 25°C  
**Detector:** UV @235nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## HINDERED PHENOLIC ANTIOXIDANTS

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** ACN/H<sub>2</sub>O 78/22 to 100/0 in 30 Min. (Linear)  
**Flow Rate:** 2.0mL/min.  
**Injection:** 100µL  
**Temperature:** 80°C  
**Detector:** UV @214nm



1. Hydroquinone monomethyl ether
2. Hydroquinone monobenzyl ether
3. 2, 5-Di-t-butylhydroquinone
4. Thiobis-(di-sec-amyphenol)
5. Tris (2-methyl-4-hydroxy-5-t-butyl phenol butane)
6. Thiobisphenol
7. 2, 2'-methylene bis (6-t-butyl-4-methyl phenol)
8. 2, 2'-methylene bis (6-t-butyl-4-ethyl phenol)
9. 4, 4'-methylene bis (2, 6-di-t-butyl phenol)
10. Pentaerythrityl tetrakis (3, 5-di-t-butyl-4-hydroxy cinnamate)
11. Octadecyl (3, 5-di-t-butyl-4-hydroxy cinnamate)




MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ANILINE, DIMETHYLANILINE AND DIETHYLANILINE

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 78/22 ACN/H<sub>2</sub>O  
**Flow Rate:** 2.0mL/min.  
**Injection:** 100µL  
**Temperature:** 25°C  
**Detector:** UV @254nm 2.0 AUFS






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## CATECHOLAMINES

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 75/24/1 0.2M NaOH/ACN/Butylamine  
**Flow Rate:** 2.0mL/min.  
**Injection:** 20µL  
**Temperature:** 25°C  
**Detector:** UV @280nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PHENOLS

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent A:** 45/45/10 H<sub>2</sub>O/ACN/MeOH w/0.1% TFA  
**Solvent B:** 8/73/19 H<sub>2</sub>O/ACN/MeOH w/0.1% TFA  
**Flow Rate:** 2.0mL/min.  
**Injection:** 5µL of 2mg/mL solution  
**Temperature:** 25°C  
**Detector:** UV @280nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ANTI-HISTAMINES

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent A:** 70/20/10 H<sub>2</sub>O/ACN/MeOH w/0.1% TFA  
**Solvent B:** 8/73/19 H<sub>2</sub>O/ACN/MeOH w/0.1% TFA  
**Flow Rate:** 2.0mL/min.  
**Injection:** 50µL of 0.13mg/mL Solution  
**Temperature:** 25°C  
**Detector:** UV @235nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ALKALOID DRUGS

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent A:** 70/20/10 H<sub>2</sub>O/ACN/MeOH w/0.1% TFA  
**Solvent B:** 10/80/10 H<sub>2</sub>O/ACN/MeOH w/0.1% TFA  
**Flow Rate:** 2.0mL/min.  
**Injection:** 100µL  
**Temperature:** 25°C  
**Detector:** UV @230nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## NUCLEOTIDES

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 97/1/2 0.01M NaOAC/ACN/MeOH  
**Flow Rate:** 2.0mL/min.  
**Injection:** 15µL  
**Temperature:** 25°C  
**Detector:** UV @254nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## NUCLEOSIDES and BASES

**Part Number:** 16010  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 10mm ID  
**Solvent:** 97/1/2 0.01M NaOAC/ACN/MeOH  
**Flow Rate:** 2.0mL/min.  
**Injection:** 50µL  
**Temperature:** 25°C  
**Detector:** UV @254nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## FATTY ACIDS (Odd Carbon)

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 95/5 MeOH/THF  
**Flow Rate:** 2.0mL/min.  
**Injection:** 100µL  
**Temperature:** 25°C  
**Detector:** UV @210nm, 0.36 AUFS








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## STEROIDS

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 85/15 ACN/H<sub>2</sub>O  
**Flow Rate:** 2.0mL/min.  
**Injection:** 50µL  
**Temperature:** 25°C  
**Detector:** UV @230nm, 0.37 AUFS



- 1. Estriol
- 2. 16-Epiestriol
- 3. 17 α-Estradiol
- 4. Esterone
- 5. Progesterone




MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## FATTY ACIDS (Even Carbon)

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 95/5 MeOH/THF  
**Flow Rate:** 2.0mL/min.  
**Injection:** 100µL  
**Temperature:** 25°C  
**Detector:** UV @210nm, 0.52 AUFS






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## SHORT CHAIN FATTY ACIDS

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 50/50 /1 ACN/H<sub>2</sub>O/H<sub>3</sub>PO<sub>4</sub>  
**Flow Rate:** 2.0mL/min.  
**Injection:** 250µL  
**Temperature:** 25°C  
**Detector:** UV @210nm, 0.6 AUFS








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## TETRACYCLINE and RELATED COMPOUNDS

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 30/70 ACN/H<sub>2</sub>O w/0.1% TFA  
**Flow Rate:** 1.5mL/min.  
**Injection:** 50µL  
**Temperature:** 25°C  
**Detector:** UV @300nm, 0.45 AUFS






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## 2,4,6-TRIS (DIMETHYLAMINOETHYL) PHENOL STANDARD (Chem Service)

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 80/20 ACN/H<sub>2</sub>O w/0.02M NaOH  
**Flow Rate:** 1.0mL/min.  
**Injection:** 5µL  
**Temperature:** 25°C  
**Detector:** UV @220nm, 1.16 AUFS






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## DIMETHYLAMINOMETHYL PHENOL STANDARD (Chem Service)

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 80/20 ACN/0.2M NaOH  
**Flow Rate:** 1.0mL/min.  
**Injection:** 10µL  
**Temperature:** 25°C  
**Detector:** UV @220nm 2.0 AUFS








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

TANNIC ACID  
(Aldrich Chemical Co)

**Part Number:** 16002  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 10cm X 10mm ID  
**Solvent:** 80/10/10 H<sub>2</sub>O/ACN/MeOH  
pH=2.5 w/H<sub>3</sub>PO<sub>4</sub>  
**Flow Rate:** 3.0mL/min.  
**Injection:** 40µL  
**Temperature:** 25°C  
**Detector:** UV @210nm 0.2 AUFS






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## GLYCEROL MONOESTERS

**Part Number:** 16507  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 78/22 to 100/0 ACN/H<sub>2</sub>O, 30 min  
**Flow Rate:** 1.0mL/min.  
**Injection:** 30µL of 1.0mg/mL  
**Temperature:** 80°C  
**Detector:** Evaporative Light Scattering






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## GLYCEROL MONOESTERS

**Part Number:** 16507  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 78/22 to 100/0 ACN/H<sub>2</sub>O, 30min  
**Flow Rate:** 1.0mL/min.  
**Injection:** 30µL of 1.0mg/mL  
**Temperature:** 80°C  
**Detector:** UV @200nm 0.2AUFS








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## SLIP AGENTS

**Part Number:** 16507  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 78/22 to 100/0 ACN/H<sub>2</sub>O 30 min.  
**Flow Rate:** 1.0mL/min.  
**Injection:** 50µL of 1.0, 0.5, 0.25 mg/mL  
**Temperature:** 80°C  
**Detector:** Evaporative Light Scattering,






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## SLIP AGENTS

**Part Number:** 16507  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 78/22 to 100/0 ACN/H<sub>2</sub>O, 30 min.  
**Flow Rate:** 1.0mL/min.  
**Injection:** 50µL of 1.0, 0.5, 0.25 mg/mL  
**Temperature:** 80°C  
**Detector:** UV @195nm, 1.0 AUFS







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

FUMONISIN B<sub>1</sub>  
(From Fusarium Moniliforme)

**Part Number:** 16508  
**Packing:** Jordi DVB Reverse Phase 1000Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 50/25/25 H<sub>2</sub>O/ACN/MeOH  
 pH=2.7 w/Formic Acid  
**Flow Rate:** 1.0mL/min.  
**Injection:** 100µL of 250ppm Solution in  
 MeOH/H<sub>2</sub>O (80/20)  
**Temperature:** 40°C  
**Detector:** Evaporative Light Scattering,








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## REVERSE PHASE STANDARD MIX (New Mobile Phase)

**Part Number:** 16507  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 50/25/25 ACN/MeOH/H<sub>2</sub>O  
**Flow Rate:** 1.0mL/min.  
**Injection:** 5µL  
**Temperature:** 50°C  
**Detector:** UV @254nm






Component	Conc (mg/mL)
1. Uracil	0.5
2. Methylparaben	0.5
3. Acetophenone	0.5
4. Butylparaben	0.5
5. Anisole	20
6. Toluene	50

Solvent: 50/25/25 ACN/MeOH/H<sub>2</sub>O  
 Temp: 50°C

Peak	Retention Time (min)	Capacity Factor	Symmetry	Efficiency (Plates/m)
1	2.91	0.00	1.4	37379
2	4.97	0.71	1.1	39051
3	7.86	1.70	1.0	40723
4	8.4	1.89	0.9	38454
5	11.33	2.90	2.7	16730
6	13.98	3.81	4.8	6857

Note: The importance of solvent/column interaction using Jordi DVB columns *cannot* be overemphasized. We have found that a 50/50/ mix of MeOH/ACN for the strong solvent is adequate for many Reverse-Phase separations and is better than either alone. We have now observed that the use of THF/ACN as strong solvent is often better than MeOH/ACN. In general Lewis bases (electron donor solvents) deactivate the aromatic rings and often dramatically increase the column efficiencies. The following data describes some Reverse Phase results for our 500Å DVB Column.






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## REVERSE PHASE STANDARD MIX (New Mobile Phase)

**Part Number:** 16507  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 30/30/40 THF/ACN/H<sub>2</sub>O  
**Flow Rate:** 1.0mL/min.  
**Injection:** 5µL  
**Temperature:** Ambient  
**Detector:** UV @254nm






Component	Conc (mg/mL)
1. Uracil	0.5
2. Methylparaben	0.5
3. Acetophenone	0.5
4. Butylparaben	0.5
5. Anisole	20
6. Toluene	50

Solvent: 30/30/40 THF/ACN/H<sub>2</sub>O  
 Temp: Ambient

Peak	Retention Time (min)	Capacity Factor	Symmetry	Efficiency (Plates/m)
1	2.68	0.00	1.6	31821
2	4.94	0.84	1.0	37592
3	6.57	1.45	0.9	40653
4	7.70	1.87	0.8	37050
5	9.81	2.65	1.8	29613
6	12.69	3.73	2.8	17307

Note: The importance of solvent/column interaction using Jordi DVB columns *cannot* be overemphasized. We have found that a 50/50/ mix of MeOH/ACN for the strong solvent is adequate for many Reverse-Phase separations and is better than either alone. We have now observed that the use of THF/ACN as strong solvent is often better than MeOH/ACN. In general Lewis bases (electron donor solvents) deactivate the aromatic rings and often dramatically increase the column efficiencies. The following data describes some Reverse Phase results for our 500Å DVB Column.









MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.


# HPLC APPLICATION

PARABENS by GPC and LC


**Part Number:** 15001  
**Packing:** Jordi DVB 500Å  
**Column:** 2-50cm X 10mm ID  
**Solvent:** See Curve detail  
**Flow Rate:** 1.5mL/min.  
**Concentration:** N/A  
**Injection:** N/A  
**Temperature:** 25°C  
**Detector:** N/A




### Strongest Deactivator



### Intermediate Deactivator



### Weakest Deactivator



Note: The importance of solvent/column interaction using Jordi DVB columns *cannot* be overemphasized. We have found that a 50/50/ mix of MeOH/ACN for the strong solvent is adequate for many Reverse-Phase separations and is better than either alone. We have now observed that the use of THF/ACN as strong solvent is often better than MeOH/ACN. In general Lewis bases (electron donor solvents) deactivate the aromatic rings and often dramatically increase the column efficiencies. The following data describes some Reverse Phase results for our 500Å DVB Column.







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## QUERCITIN and RUTIN

**Part Number:** 16011  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 10mm ID  
**Solvent:** 75/25 Acetic Acid/MeOH  
**Flow Rate:** 2.0mL/min.  
**Injection:** 100µL  
**Temperature:** Ambient  
**Detector:** UV 360nm, 2 AUFS






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## QUERCITIN and RUTIN

**Part Number:** 15001  
**Packing:** Jordi DVB 500Å  
**Column:** 50cm X 10mm ID  
**Mobile Phase:** 50/50 Pyridine/Glacial Acetic Acid  
**Flow Rate:** 1.0mL/min.  
**Injection:** 15µL  
**Temperature:** 40°C  
**Detector:** UV @ 365nm, Sens. 0.2AUFS,  
 Press. 2500 PSIG

Component	Conc (mg/mL)
1. Rutin	0.1
2. Quercitin	0.1



Peak	Retention Time (min)	Capacity Factor	Symmetry	Efficiency (Plates/m)
1	26.11	0.00	1.4	46816
2	29.88	0.14	1.3	54569





MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## 2'DEOXY-2'FLUOROADENOSINE ALPHATHIOPHOSPHATE

**Part Number:** 16507  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 4.6mm ID  
**Gradient:** The elution order of the standards is the monophosphate, the diastereomer triphosphate and finally the free nucleoside  
**Buffer A:** 0.1M TEAA pH 7.0/ Acetonitrile 99/1  
**Buffer B:** 0.1M TEAA pH 7.0/ Acetonitrile/ Methanol 50/25/25  
**Flow Rate:** 1.0mL/min.  
**Injection:** 10µL 10mM Standards  
**Temperature:** Ambient  
**Detector:** UV @260nm



	Gradient Linear				
Time (min)	0	10	11	12	15
%A	90	0	0	90	90
%B	10	100	100	10	10

This column readily resolves this  $\alpha$ -thiophosphate diastereomer. By switching to methanol versus a mixture of acetonitrile and methanol, higher baseline resolution can be achieved, even between the diastereomers.







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## 2'DEOXY-2'FLUOROADENOSINE $\alpha$ -THIOTRIPHOSPHATE, SODIUM SALT

**Part Number:** 16507  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 4.6mm ID  
**Gradient:** The elution order of the standards is the monophosphate, the diastereomer triphosphate and finally the free nucleoside  
**Buffer A:** 0.1M TEAA pH 7.0/ Acetonitrile 99/1  
**Buffer B:** 0.1M TEAA pH 7.0/ Acetonitrile/ Methanol 50/25/25  
**Flow Rate:** 1.0mL/min.  
**Injection:** 10 $\mu$ L 10mM Standards  
**Temperature:** Ambient  
**Detector:** UV @260nm



Time (min)	Gradient Linear				
	0	10	11	12	15
%A	90	0	0	90	90
%B	10	100	100	10	10

This column readily resolves this  $\alpha$ -thiophosphate diastereomer. By switching to methanol versus a mixture of acetonitrile and methanol, higher baseline resolution can be achieved, even between the diastereomers.




MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## 2'DEOXY-2'FLUOROADENOSINE

**Part Number:** 16507  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 4.6mm ID  
**Gradient:** The elution order of the standards is the monophosphate, the diastereomer triphosphate and finally the free nucleoside  
**Buffer A:** 0.1M TEAA pH 7.0/ Acetonitrile 99/1  
**Buffer B:** 0.1M TEAA pH 7.0/ Acetonitrile/ Methanol 50/25/25  
**Flow Rate:** 1.0mL/min.  
**Injection:** 10µL 10mM Standards  
**Temperature:** Ambient  
**Detector:** UV @260nm



	Gradient Linear				
Time (min)	0	10	11	12	15
%A	90	0	0	90	90
%B	10	100	100	10	10

This column readily resolves this  $\alpha$ -thiophosphate diastereomer. By switching to methanol versus a mixture of acetonitrile and methanol, higher baseline resolution can be achieved, even between the diastereomers.




MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## MIX of 2'DEOXY-2'FLUOROADENOSINE ALPHATHIOTRIPHOSPHATE and 2'DEOXY-2'FLUOROADENOSINE

**Part Number:** 16507  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 4.6mm ID  
**Gradient:** The elution order of the standards is the monophosphate, the diastereomer triphosphate and finally the free nucleoside  
**Buffer A:** 0.1M TEAA pH 7.0/ Acetonitrile 99/1  
**Buffer B:** 0.1M TEAA pH 7.0/ Acetonitrile/  
Methanol 50/25/25  
**Flow Rate:** 1.0mL/min.  
**Injection:** 10µL  
**Temperature:** Ambient  
**Detector:** UV @260nm



This column readily resolves this  $\alpha$ -thiophosphate diastereomer. By switching to methanol versus a mixture of acetonitrile and methanol, higher baseline resolution can be achieved, even between the diastereomers.






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

MIX of 2'DEOXY-2'FLUOROADENOSINE and the MONO and TRI-

**Part Number:** 16507  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 4.6mm ID  
**Gradient:** The elution order of the standards is the monophosphate, the diastereomer triphosphate and finally the free nucleoside  
**Buffer A:** 0.1M TEAA pH 7.0/ Acetonitrile 99/1  
**Buffer B:** 0.1M TEAA pH 7.0/ Acetonitrile/ Methanol 50/25/25  
**Flow Rate:** 1.0mL/min.  
**Injection:** 10µL  
**Temperature:** Ambient  
**Detector:** UV @260nm



This column readily resolves this  $\alpha$ -thiophosphate diastereomer. By switching to methanol versus a mixture of acetonitrile and methanol, higher baseline resolution can be achieved, even between the diastereomers.




MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## MIX of 2'DEOXY-2'FLUOROADENOSINE ALPHATHIOTRIPHOSPHATE and 2'DEOXY-2'FLUOROADENOSINE

**Part Number:** 16507  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 25cm X 4.6mm ID  
**Buffer A:** 99/1 0.1M TEAA, pH 7.0/Acetonitrile  
**Buffer B:** 50/50 0.1M TEAA, pH 7.0/Acetonitrile  
**Flow Rate:** 1.0mL/min.  
**Injection:** 10µL  
**Temperature:** Ambient  
**Detector:** UV @260nm



This column readily resolves this  $\alpha$ -thiophosphate diastereomer. By switching to methanol versus a mixture of acetonitrile and methanol, higher baseline resolution can be achieved, even between the diastereomers.




MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

RP MIX

**Part Number:** 60026  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 53cm X 7mm ID  
**Solvent:** 10/20/30/40 MeOH/ACN/THF/H<sub>2</sub>O  
**Flow Rate:** 1.0mL/min.  
**Injection:** 5µL  
**Concentration:** 0.5mg/mL  
**Temperature:** 50°C  
**Detector:** UV @254nm, Sens. 2.0 AUFS, Press. 650 psig








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

SEPARATION of POLYETHYLENE GLYCOL and POLYPROPYLENE GLYCOL

**Part Number:** 16502  
**Packing:** Jordi DVB Reverse Phase 500Å  
**Column:** 3- 15cm X 4.6mm ID  
**Gradient:** 40/60 A/B → 100 B linear over 10 min.  
A: 78/22 CH<sub>3</sub>CN/H<sub>2</sub>O B: 2-propanol  
**Flow Rate:** 1.0mL/min.  
**Injection:** 20µL  
**Concentration:** 1mg/mL  
**Temperature:** 80°C  
**Detector:** Alltech Mark III ELSD, Attn: 16, Exhaust: 50°C, Flow 2.00 SLPM Air









MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## SEPARATION of POLYETHYLENE GLYCOL and POLYPROPYLENE GLYCOL

**Part Number:** 16502  
**Packing:** Jordi DVB 500Å  
**Column:** 3- 15cm X 4.6mm ID  
**Isocratic:** 2-propanol  
**Flow Rate:** 1.0mL/min.  
**Injection:** 20µL  
**Concentration:** 1mg/mL  
**Temperature:** 80°C  
**Detector:** Alltech Mark III ELSD, Attn: 16,  
 Exhaust: 50°C, Flow 2.00 SLPM Air







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ANTIULCERATIVE OMEPRAZOLE

**Part Number:** 16502  
**Packing:** Jordi DVB RP 500Å  
**Column:** 15cm X 4.6mm ID  
**Mobile Phase:** 10/20/30/40 MeOH/ACN/THF/H<sub>2</sub>O  
**Flow Rate:** 1.0mL/min.  
**Injection:** 20µL  
**Concentration:** 1mg/mL  
**Temperature:** 50°C  
**Detector:** UV @254nm, 2.0AUFS








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

PEG 425

**Part Number:** 16508  
**Packing:** Jordi DVB Reverse Phase 10<sup>3</sup>Å  
**Column:** 25cm x 4.6mm  
**Solvent:** 60/40MeOH/ H<sub>2</sub>O  
**Flow Rate:** 0.5mL/min.  
**Injection:** 50µL  
**Temperature:** N/A  
**Detector:** Alltech ELSD Attn 16, SLPM 16, Drift Tube 80°C






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## DILAURYLTHIODIPROPIONATE, DIMYRISTYLTHIODIPROPIONATE, & DISTEARYLTHIODIPROPIONATE-POLYMER ANTIOXIDANTS

**Part Number:** 15310  
**Packing:** Jordi DVB 10<sup>3</sup>Å  
**Column:** 25cm X 4.6mm ID  
**Gradient:** 100% ACN → 100% 45/45/10 ACN/IPA/IO linear over 30 min  
**Flow Rate:** 1.0mL/min.  
**Injection:** 10µL (1mg/mL ea. std. in 50/50 CHCl<sub>3</sub>/IO)  
**Temperature:** 80°C  
**Detector:** Waters 990 PDA & Alltech 500 ELSD  
**990 PDA Cond:** 210nm 1 AUFS  
**ELSD Conditions:** Attn 16, Drift Tube 80°C, 2.0 SLPM House Air, Time 1 sec.






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## IRGANOX 1010, NAUGARD 445 & DISTEARYLTHIODIPROPIONATED-POLYMER ANTIOXIDANTS

**Part Number:** 15310  
**Packing:** Jordi DVB 10<sup>3</sup>Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 100% ACN → 100% 45/45/10 ACN/IPA/IO linear over 30min.  
**Flow Rate:** 1.0mL/min.  
**Injection:** 10µL (1mg/mL ea. std. in 50/50 CHCl<sub>3</sub>/IO)  
**Temperature:** 80°C  
**Detector:** Waters 990 PDA & Alltech 500 ELSD  
**Waters 990 Cond.:** 215nm 1 AUFS  
**Alltech ELSD Cond:** Attn 16, Drift Tube 80°C, 2.0 SLPM House Air, Time 1 sec.








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## REVERSE PHASE MIX

**Part Number:** 16508  
**Packing:** Jordi DVB Reverse Phase 10<sup>3</sup>Å  
**Column:** 25cm X 4.6mm ID  
**Mobile Phase:** 10/20/30/40 MeOH/ACN/THF/H<sub>2</sub>O  
**Flow Rate:** 1.0mL/min.  
**Injection:** 3µL  
**Conc.:** 0.5mg/mL of Uracil, Methylparaben, Propylparaben, Butylparaben  
**Temperature:** 50°C  
**Detector:** UV @254nm, Sens. 2.0 AUFS






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.


# HPLC APPLICATION

## SEA NINE 211 ROHM & HAAS ANTIFOULING AGENT

**Part Number:** 15310  
**Packing:** Jordi DVB Organic 10<sup>3</sup>Å  
**Column:** 25cm X 4.6mm ID  
**Gradient:** 78/22 → 100/0 CH<sub>3</sub>CN/H<sub>2</sub>O linear over 30min.  
**Flow Rate:** 1.0mL/min.  
**Injection:** See Curve for details  
**Temperature:** 80°C  
**Detector:** Waters 990 PDA, 225nm or 283nm 1.0 AUFS



10 $\mu$ L inj. 1/100 dilution of Sea Nine 211-30% (4,5-dichloro-2-n-octyl-4-isothiazolin-3-one and 70% xylene). Thus each injection contained approximately 7mg xylene and 3mg active ingredient






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## URETHANE PREPOLYMER

**Part Number:** 16508  
**Packing:** Jordi DVB Reverse Phase 10<sup>3</sup>Å  
**Column:** 25cm X 4.6mm ID  
**Gradient:** 100/0 → 0/100 ACN/CHCl<sub>3</sub> linear over 30min.  
**Flow Rate:** 1.0mL/min.  
**Injection:** 5μL  
**Temperature:** 65°C  
**Detector:** Alltech 500 ELSD, see curve for details  
**Conc.:** 100mg/mL









MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## HYDROQUINONE

**Part Number:** 16508  
**Packing:** Jordi DVB Reverse Phase 10<sup>3</sup>Å  
**Column:** 25cm X 4.6mm ID  
**Isocratic:** 10/10/80 MeOH/ACN/H<sub>2</sub>O w/0.1% TFA for 15min. then flush with THF for 15min. and then re-equilibrate with initial mobile phase.  
**Flow Rate:** 1.0mL/min.  
**Injection:** 100µL  
**Temperature:** 80°C  
**Detector:** Waters 990 PDA @290nm  
**Conc.:** 10µg/mL






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## POLYMER SAMPLE SEPARATIONS USING RP CHROMATOGRAPHY

**Part Number:** 16508  
**Packing:** Jordi DVB Reverse Phase 10<sup>3</sup>Å  
**Column:** 25cm X 4.6mm ID  
**Gradient:** 100/0 → 0/100 ACN/CHCl<sub>3</sub> linear over 30 min.  
**Flow Rate:** 1.0mL/min.  
**Injection:** 20µL  
**Conc.:** 5mg/mL  
**Temperature:** 65°C  
**Detector:** Alltech 500 ELSD: ATTN 16, Heater 65°C  
Exhaust 43°C, Flow 2.00 SLPM Air






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ANTI STATIC ANALYSIS

**Part Number:** 16508  
**Packing:** Jordi DVB Reverse Phase $10^3\text{\AA}$   
**Column:** 25cm X 4.6mm ID  
**Gradient:** 78/22  $\rightarrow$  100/0  $\text{CH}_3\text{CN}/\text{H}_2\text{O}$  over 30 min. linear  
**Flow Rate:** 1.0mL/min.  
**Injection:** 10 $\mu\text{L}$   
**Concentration:** 4mg/mL  
**Temperature:** 80 $^\circ\text{C}$   
**Detector:** Waters 990 PDA, see curve for conditions










MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## MINERAL and MOTOR OILS

**Part Number:** 16508  
**Packing:** Jordi DVB RP 10<sup>3</sup>Å  
**Column:** 25cm X 4.6mm ID  
**Isocratic:** 75/25 Ethanol/2-Propanol  
**Flow Rate:** 1.0mL/min.  
**Injection:** See Curve Detail  
**Temperature:** 80°C  
**Detector:** Alltech 500 ELSD  
**ELSD Conditions:** ATTN 4 Heater 70°C Exhaust 50°C Flow 2.0 SLPM






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION



## NYLON 6 MONOMER ANALYSIS

**Part Number:** 16508  
**Packing:** Jordi DVB Reverse Phase10<sup>3</sup>Å  
**Column:** 25cm X 4.6mm ID  
**Gradient:** 10/90 → 100/0 CH<sub>3</sub>CN/H<sub>2</sub>O  
 Linear over 30 min.  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40µL of 10mg/mL extract  
**Temperature:** 80°C  
**Detector:** Waters 990 PDA, see  
 curve detail for conditions



NYLON SPOON EXTRACT

**ε-caprolactam**






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## BENZOIC ACID

**Part Number:** 16508  
**Packing:** Jordi Gel DVB  
**Column:** 1-250 X 4.6 mm 10<sup>3</sup> Å  
**Solvent:** 35/5/60 MeOH/THF/H<sub>2</sub>O+0.1% TFA  
**Flow Rate:** 1.0mL/min.  
**Injection:** 10µL  
**Concentration:** 1 mg/mL Standards 10 µg/inj each Standard  
**Temperature:** 80°C  
**Detector:** Waters 990 PDA , 240 nm 0.6 AuFS








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PROTEINS and ENZYMES

**Part Number:** 17020  
**Packing:** Jordi DVB Solid Bead  
**Column:** 10cm X 10mm ID  
**Solvent A:** 0.1% TFA in H<sub>2</sub>O  
**Solvent B:** 0.1% TFA in ACN  
**Flow Rate:** 3.0mL/min.  
**Injection:** 20µL  
**Temperature:** 25°C  
**Detector:** UV @215nm






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PROTEINS

**Part Number:** 17020  
**Packing:** Jordi DVB Solid Bead  
**Column:** 10cm X 10mm ID  
**Solvent A:** 0.1% TFA in H<sub>2</sub>O  
**Solvent B:** 0.1% TFA in ACN  
**Flow Rate:** 3.0mL/min.  
**Injection:** 50µL  
**Temperature:** 25°C  
**Detector:** UV @215nm





# Glucose DVB








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PHOSPHOLIPID MIXTURE

**Part Number:** 15061  
**Packing:** Jordi DVB Glucose 500Å  
**Column:** 25cm X 10mm ID  
**Solvent:** 48/42/10 CHCl<sub>3</sub>/MeOH/H<sub>2</sub>O w/0.2% NH<sub>4</sub>OH  
**Flow Rate:** 1.0mL/min.  
**Injection:** 2µL of 2mg/mL solution  
**Temperature:** 25°C  
**Detector:** Evaporative Light Scattering, Varex IIA






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## LIPIDS and A TRIGLYCERIDE

**Part Number:** 18610  
**Packing:** Jordi DVB Glucose NP 500Å  
**Column:** 25cm X 10mm ID  
**Mobile Phase:** 50/43/7 Chloroform/Methanol/0.15% TFA in H<sub>2</sub>O  
**Flow Rate:** 1.0mL/min.  
**Injection:** N/A  
**Concentration:** N/A  
**Temperature:** Ambient  
**Detector:** ELSD MKIII, Drift Tube Temp. 80°C with Nitrogen @ 2.20 SLPM





# Hydroxylate DVB









MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## ANTIULCERATIVE OMEPRAZOLE

**Part Number:** 21000  
**Packing:** Jordi DVB Hydroxylated RP 500Å  
**Column:** 15cm X 4.6mm ID  
**Mobile Phase:** 10/20/20/50 MeOH/ACN/THF/H<sub>2</sub>O  
**Flow Rate:** 1.0mL/min.  
**Injection:** 5.0µL  
**Concentration:** 1mg/mL  
**Temperature:** 50°C  
**Detector:** UV @254nm, 1.0AUFS



Peak	Retention Time(min)	Capacity Factor	Symmetry	Efficiency (Plates/m)
1	2.92	0.00	1.8	16232



# Organic Acid DVB








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## SEPARATION of ORGANIC ACIDS

**Part Number:** 17000  
**Packing:** Jordi DVB Organic Acid 500Å  
**Column:** 50cm X 10mm ID  
**Solvent :** 0.01M Phosphoric Acid, pH 3 w/NaOH  
**Flow Rate:** 1.5mL/min.  
**Injection:** 200µL  
**Concentration:** 10mg/mL  
**Temperature:** Ambient  
**Detector:** UV -210nm, RI @4X





# Polyamino DVB








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## SUGAR STANDARDS

**Part Number:** 17010  
**Packing:** Jordi DVB Polyamino 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 80/15/5 ACN/H<sub>2</sub>O/MeOH  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40µL  
**Concentration:** 2mg/mL  
**Temperature:** Ambient  
**Detector:** RI






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.



# HPLC APPLICATION

## AZODICARBONAMIDE

**Part Number:** 60055  
**Packing:** Jordi DVB Polyamino RP 500Å  
**Column:** 5.3cm X 7mm ID  
**Solvent:** 90/10 H<sub>2</sub>O/DMSO  
**Flow Rate:** 1.0mL/min.  
**Injection:** 10µL  
**Concentration:** 1mg/mL  
**Temperature:** Ambient  
**Detector:** Waters 990 PDA, 265nm



ADCA








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## POLYHYDROXY ALKALOIDS

**Part Number:** 17051  
**Packing:** Jordi NP Polyamino 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent A:** Acetonitrile:0.2%(w/v)  
**Solvent B:** Ammonium Acetate in H<sub>2</sub>O(15:85)  
**Flow Rate:** 0.5mL/min.  
**Injection:** 5µg HNJ  
6.25µg DMDP  
**Temperature:** Ambient  
**Detector:** ELSD MKIII






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## SUGAR STANDARDS

**Part Number:** 17010  
**Packing:** Jordi DVB Polyamino 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 75/20/5 ACN/H<sub>2</sub>O/MeOH  
**Flow Rate:** 1.0mL/min.  
**Injection:** 10µL of 2% Concentration  
**Det. Temperature:** 30°C  
**Detector:** RI,4X






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

CORN SYRUP by Gradient Elution

**Part Number:** 17010  
**Packing:** Jordi DVB Polyamino 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent A:** 85/15 ACN/H<sub>2</sub>O  
**Solvent B:** H<sub>2</sub>O  
**Gradient:** 15 min Linear, 0%B to 21%B  
**Flow Rate:** 1.0mL/min.  
**Injection:** 100µL of 1% Concentration  
**Temperature:** 25°C  
**Detector:** Evaporative Light Scattering, Varex ELSD II @ 100X





# Peptide Protein DVB






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PROTEIN STANDARD SEPARATION

**Part Number:** 10001  
**Packing:** Jordi Peptide/Protein 10<sup>3</sup>Å  
**Column:** 5cm X 4.6mm ID  
**Gradient:** 80/20 → 40/60 A/B over 30 min. linear  
**Solvent A:** 0.15% TFA in H<sub>2</sub>O  
**Solvent B:** 0.15% TFA in ACN  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40µL  
**Concentration:** Protein Stds. diluted to 4mL in 75/25 A/B  
**Temperature:** Ambient  
**Detector:** Waters 990, 220nm








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## JORDI PEPTIDE PROTEIN

### Polyamide Co-Polymer 5 Peptide Standards

**Part Number:** 60060  
**Packing:** Jordi Peptide Protein 10<sup>4</sup>Å  
**Column:** 5.3cm X 7mm ID  
**Solvent:** 90/10 → 60/40 A/B over 30 min. Linear  
A: 99/1 H<sub>2</sub>O/Acetic Acid B: 99/1 ACN/Acetic Acid  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40µL  
**Concentration:** HPLC Pesticide Standards diluted to 4mL in 99/1/H<sub>2</sub>O/Acetic Acid  
**Temperature:** Ambient  
**Detector:** UV-270nm @0.025 AUFS








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PROTEIN STANDARD SEPARATION

**Part Number:** 10001  
**Packing:** Jordi Peptide Protein 10<sup>4</sup>Å  
**Column:** 5cm X 4.6mm ID  
**Gradient:** 90/10 → 60/40 A/B over 30 min. linear  
**Solvent A:** 98/2 H<sub>2</sub>O/Acetic Acid  
**Solvent B:** 98/2 ACN/Acetic Acid  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40µL  
**Concentration:** Protein Standards diluted to 4mL in 75/25 A/B  
**Temperature:** Ambient  
**Detector:** UV-270nm @ 0.01 AUFS







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PROTEIN STANDARD SEPARATION

**Part Number:** 10001  
**Packing:** Jordi Peptide Protein 10<sup>4</sup>Å  
**Column:** 5cm X 4.6mm ID  
**Gradient:** 80/20 → 40/60 A/B over 10 min. linear  
**Solvent A:** 0.15% TFA in H<sub>2</sub>O  
**Solvent B:** 0.15% TFA in ACN  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40µL  
**Concentration:** Protein Standards diluted to 4mL in 75/25 A/B  
**Temperature:** Ambient  
**Detector:** UV (see curve for details)






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PROTEIN STANDARD SEPARATION

**Part Number:** 10001  
**Packing:** Jordi Peptide Protein 10<sup>4</sup>Å  
**Column:** 5cm X 4.6mm ID  
**Gradient:** 80/20 → 40/60 A/B over 10 min. linear  
**Solvent A:** 0.15% TFA in H<sub>2</sub>O  
**Solvent B:** 0.15% TFA in ACN  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40µL  
**Concentration:** Protein Standards diluted to 4mL in 75/25 A/B  
**Temperature:** Ambient  
**Detector:** UV-270nm @0.03AUFS











MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PROTEIN STANDARD SEPARATION

**Part Number:** 10001  
**Packing:** Jordi Peptide Protein 10<sup>4</sup>Å  
**Column:** 5cm X 4.6mm ID  
**Gradient:** 80/20 → 40/60 A/B over 10 min. linear  
**Solvent A:** 0.15% TFA in H<sub>2</sub>O  
**Solvent B:** 0.15% TFA in ACN  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40µL  
**Concentration:** 0.1mg/mL  
**Temperature:** Ambient  
**Detector:** UV-220nm @0.3AUFS









MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PROTEIN STANDARD SEPARATION

**Part Number:** 10001  
**Packing:** Jordi Peptide Protein 10<sup>4</sup>Å  
**Column:** 5cm X 4.6mm ID  
**Gradient:** 80/20 → 40/60 A/B over 10 min. linear  
**Solvent A:** 0.15% TFA in H<sub>2</sub>O  
**Solvent B:** 0.15% TFA in ACN  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40µL  
**Concentration:** 0.1mg/mL  
**Temperature:** Ambient  
**Detector:** UV-220nm@ 0.3AUFS









MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PROTEIN STANDARD SEPARATION

**Part Number:** 10001  
**Packing:** Jordi Peptide Protein 10<sup>4</sup>Å  
**Column:** 5cm X 4.6mm ID  
**Gradient:** 80/20 → 40/60 A/B over 10 min. linear  
**Solvent A:** 0.15% TFA in H<sub>2</sub>O  
**Solvent B:** 0.15% TFA in ACN  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40µL  
**Concentration:** 0.1mg/mL  
**Temperature:** Ambient  
**Detector:** UV-220nm@ 0.3AUFS











MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PROTEIN STANDARD SEPARATION

**Part Number:** 10001  
**Packing:** Jordi Peptide Protein 10<sup>4</sup>Å  
**Column:** 5cm X 4.6mm ID  
**Gradient:** 80/20 → 40/60 A/B over 10 min. linear  
**Solvent A:** 0.15% TFA in H<sub>2</sub>O  
**Solvent B:** 0.15% TFA in ACN  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40µL  
**Concentration:** 0.1mg/mL  
**Temperature:** Ambient  
**Detector:** UV-220nm@ 0.3AUFS






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PEPTIDE PROTEIN STANDARD SEPARATION

**Part Number:** 10001  
**Packing:** Jordi Peptide Protein 10<sup>4</sup>Å  
**Column:** 2-5cm X 4.6mm ID  
**Gradient:** 80/20 → 40/60 A/B over 30 min. linear  
**Solvent A:** 0.15% TFA in H<sub>2</sub>O  
**Solvent B:** 0.15% TFA in ACN  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40μL Protein Standard Mix  
40μL HPLC Peptide Standard Mix  
**Concentration:** Protein Standards diluted to 4mL in 75/25 A/B  
**Temperature:** Ambient  
**Detector:** N/A







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## BOVINE SERUM ALBUMIN

**Part Number:** 10001  
**Packing:** Jordi Peptide Protein 10<sup>4</sup>Å  
**Column:** 5cm X 4.6mm ID  
**Gradient:** 80/20 → 40/60 A/B over 30 min. linear  
**Solvent A:** 0.15% TFA in H<sub>2</sub>O  
**Solvent B:** 0.15% TFA in ACN  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40μL  
**Concentration:** Protein Standards Diluted to 4mL in 75/25 A/B @ 1mg/mL  
**Temperature:** Ambient  
**Detector:** UV-270nm @ 0.05AUFS UV-220nm @ 0.3AUFS








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PROTEIN STANDARDS MIX

**Part Number:** 10001  
**Packing:** Jordi Peptide Protein 10<sup>4</sup>Å  
**Column:** 5cm X 4.6mm ID  
**Gradient:** 80/20 → 40/60 A/B over 30 min. linear  
**Solvent A:** 0.15% TFA in H<sub>2</sub>O  
**Solvent B:** 0.15% TFA in ACN  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40µL  
**Concentration:** Protein Standards Diluted to 4mL in 75/25 A/B  
**Temperature:** Ambient  
**Detector:** UV-270nm @ 0.03AUFS  
UV-220nm @ 0.3AUFS






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PROTEIN STANDARD SEPARATION

**Part Number:** 10001  
**Packing:** Jordi Peptide Protein 10<sup>4</sup>Å  
**Column:** 5cm X 4.6mm ID  
**Gradient:** 80/20 → 40/60 A/B over 30 min. linear  
**Solvent A:** 0.15% TFA in H<sub>2</sub>O  
**Solvent B:** 0.15% TFA in ACN  
**Flow Rate:** 1.0mL/min.  
**Injection:** 40μL  
**Concentration:** Protein Standards Diluted to 4mL in 75/25 A/B  
**Temperature:** Ambient  
**Detector:** UV-220nm, 0.5 AUFS





# Sulfonated DVB








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## DIVINYLBENZENE

**Part Number:** 15252  
**Packing:** Jordi DVB Sulfonated 500Å in Ag<sup>+</sup> form  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 85/15 CHCl<sub>3</sub>/MeOH  
**Flow Rate:** 0.7mL/min.  
**Injection:** 3µL  
**Temperature:** 25°C  
**Detector:** UV @260nm, 1.6AUFS







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## PHOSPHORIC COMPOUNDS

**Part Number:** 15254  
**Packing:** Jordi DVB Sulfonated 500Å  
**Column:** 25cm X 4.6mm ID  
**Solvent:** 0.02M Phosphoric Acid w/3%ACN  
 pH 2.5w/ NaOH  
**Flow Rate:** 1.0mL/min.  
**Injection:** 50µL of 1mg/mL Solution  
**Temperature:** 25°C  
**Detector:** RI






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## GLUCOSE & ACETIC ACID

**Part Number:** 15041  
**Packing:** Jordi DVB Sulfonated 500Å  
**Column:** 25cm X 10mm ID  
**Mobile Phase:** 0.01M Phosphoric Acid Adj. to pH 7.5  
**Flow Rate:** 2.0mL/min.  
**Injection:** 20µL  
**Temperature:** Ambient  
**Detector:** RI 4X





# Sax Quat

## DVB









MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## CHROMATE PLATING BATH

**Part Number:** 18701  
**Packing:** Jordi DVB SAX Quat 10<sup>3</sup>Å  
**Column:** 10cm X 4.6mm ID  
**Solvent:** 90/10 15mL/L EZ Lute Concentrate/MeOH  
**Flow Rate:** 1.0mL/min.  
**Injection:** 20µL  
**Temperature:** 40°C  
**Detector:** Waters 430 Conductivity Gain 0.5,  
 Range 0.5 0.25µS FS







MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## CORROSION TEST

**Part Number:** 18701  
**Packing:** Jordi DVB SAX Quat 10<sup>3</sup>Å  
**Column:** 10cm X 4.6mm ID  
**Solvent:** 90/10 15mL/L EZ Lute Concentrate/MeOH  
**Flow Rate:** 1.0mL/min.  
**Injection:** 20µL  
**Temperature:** 40°C  
**Detector:** Waters 430 Conductivity Gain 0.2, Range 0.5 1µS FS








MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## 7 ANION STANDARD MIX

**Part Number:** 18701  
**Packing:** Jordi DVB SAX Quat 10<sup>3</sup>Å  
**Column:** 10cm X 4.6mm ID  
**Solvent:** 10mL/L EZ Lute  
**Flow Rate:** 9.0mL/min.  
**Injection:** 20µL  
**Temperature:** Ambient  
**Detector:** Conductivity






MATERIAL SOLUTIONS. UNCOMPROMISING INTEGRITY.

# HPLC APPLICATION

## 7 ANION STANDARD MIX

**Part Number:** 18701  
**Packing:** Jordi DVB SAX Quat 10<sup>3</sup>Å  
**Column:** 10cm X 4.6mm ID  
**Solvent:** 10mL/L EZ Lute  
**Flow Rate:** 1.0mL/min.  
**Injection:** 20µL  
**Temperature:** Ambient  
**Detector:** Conductivity



# INDEX

Page	Compound	Solvent	Column
Pg 26	Alkaloid Drugs	70/20/10 H <sub>2</sub> O/ACN/MeOH w/0.1% TFA	DVB C18
Pg 61	Alkaloid Drugs	A: 70/20/10 B: 10/80/10 H <sub>2</sub> O/ACN/MeOH w/0.1% TFA	DVB/RP
Pg 32	Ampicillin	50/50 ACN/0.05M KH <sub>2</sub> PO <sub>4</sub>	DVB C18
Pg 57	Aniline, Dimethylaniline & Diethylaniline	78/22 ACN/H <sub>2</sub> O	DVB/RP
Pg 140	7 Anion Std. Mix @9mL/min	10mL/L EZ Lute Concentrate	DVB SAX Quat
Pg 141	7 Anion Std. Mix @1mL/min	10mL/L EZ Lute Concentrate	DVB SAX Quat
Pg 17	Anticonvulsants	A: H <sub>2</sub> O B: 87.5/12.5 ACN/MeOH	DVB C18
Pg 60	Antihistamines	A: 70/20/10 B: 8/73/19 H <sub>2</sub> O/ACN/MeOH w/0.1% TFA	DVB/RP
Pg 100	Antistatic Analysis	78/22→100/0 CH <sub>3</sub> CN/H <sub>2</sub> O	DVB RP 103Å
Pg 110	Antiulcerative Omperazole	10/20/20/50 MeOH/ACN/THF/H <sub>2</sub> O	DVB Hydroxylated RP 500Å
Pg 91	Antiulcerative Omperazole	10/20/30/40 MeOH/ACN/THF/H <sub>2</sub> O	DVB RP 500Å
Pg 37	Ascorbic and Dehydroascorbic Acid	H <sub>2</sub> O, pH 2.7 w/acetic acid: Acetonitrile/Methanol (5:1) 40:60	DVB C18
Pg 53	Aspirin & Related Compounds	50/50 ACN/H <sub>2</sub> O w/0.5% TFA	DVB/RP
Pg 8	Aspirin & Related Compounds	60/40 ACN/H <sub>2</sub> O w/0.1% TFA	DVB C18
Pg 21	Aspirin & Related Compounds	40/60 ACN/H <sub>2</sub> O w/0.1% TFA	DVB C18
Pg 115	Azodicarbonamide	90/10 H <sub>2</sub> O/DMSO	DVB Polyamide RP 500Å
Pg 18	Barbiturates	70/22.5/7.5 0.05M Na <sub>2</sub> HPO <sub>4</sub> /ACN/MeOH	DVB C18
Pg 48	Base-Neutral-Acid Test Mix	55/37.2/7.8 ACN/H <sub>2</sub> O/MeOH	DVB/RP
Pg 103	Benzoic Acid	35/5/60 MeOH/THF/H <sub>2</sub> O+0.1% TFA	DVB/RP
Pg 130	Bovine Serum Albumin	80/20→40/60 A/B 30 Min. Linear A: 0.15% TFA in H <sub>2</sub> O B: 0.15% TFA in ACN	Jordi Peptide/ Protein Hydroxylated
Pg 9	Butyl Phenol Standards	70/30 ACN/H <sub>2</sub> O w/ 0.1% TFA	DVB C18
Pg 49	Butylphenol Standards	80/20 ACN/H <sub>2</sub> O w/0.05% TFA	DVB/RP
Pg 45	Caffeine in Coffee	50/15/35 0.01M LiNO <sub>3</sub> /ACN/MeOH	DVB/RP
Pg 47	Caffeine in Coke & Pepsi	50/15/35 0.01M LiNO <sub>3</sub> /ACN/MeOH	DVB/RP
Pg 46	Caffeine in Tea & Chocolate	50/15/35 0.01M LiNO <sub>3</sub> /ACN/MeOH	DVB/RP
Pg 58	Catecholamines	75/24/1 0.2M NaOH/ACN/Butylamine	DVB/RP
Pg 25	Catecholamines	A: 98/1/1 B: 79/20/1 0.01M NaH <sub>2</sub> PO <sub>4</sub> /Butylamine	DVB C18
Pg 19	Choline Compounds	H <sub>2</sub> O @pH 3.3 w/HCl	DVB C18
Pg 10	Choline Compounds	90/3/7 0.1M Na <sub>2</sub> HPO <sub>4</sub> /ACN/MeOH	DVB C18
Pg 138	Chromate Plating Bath	90/10 15mL/L EZ Lute Conc./MeOH	DVB Anion SAX Quat
Pg 54	Column Test Mix	75/20/5 ACN/H <sub>2</sub> O/MeOH w/0.1% TFA	DVB/RP
Pg 11	Column Test Mix	65/30/5 ACN/H <sub>2</sub> O/MeOH w/0.1% TFA	DVB C18
Pg 118	Corn Syrup	A: 85/15 ACN/H <sub>2</sub> O B: H <sub>2</sub> O	DVB Polyamine
Pg 139	Corrosion Test	90/10 15mL/L EZ Lute Conc./MeOH	DVB Sax Quat
Pg 84	2'-deoxy-2'-fluoroadenosine	A: 0.1M TEAA/ACN 99:1 B: 0.1M TEAA/ACN/MeOH 50:25:25	DVB/RP
Pg 86	Mix of 2'-deoxyadenosine & the Mono & Tri	A: 0.1M TEAA/ACN 99:1 B: 0.1M TEAA/ACN/MeOH 50:25:25	DVB/RP
Pg 82	2'-deoxy-2'-fluoroadenosine alphathiosphosphate	A: 0.1M TEAA/ACN 99:1 B: 0.1M TEAA/ACN/MeOH 50:25:25	DVB/RP
Pg 85	2'-deoxy-2'-fluoroadenosine alphatrithiosphosphate & 2'-deoxy-2'-fluoroadenosine	A: 0.1M TEAA/ACN 99:1 B: 0.1M TEAA/ACN/MeOH 50:25:25	DVB/RP



Pg 87	2'-deoxy-2'-fluoroadenosine alphatrithiophosphate & 2'-deoxy-2'-fluoroadenosine	A: 0.1M TEAA/1% ACN B: 50% 0.1M TEAA/ACN	DVB/RP
Pg 83	2'-deoxy-2'-fluoroadenosine alphatrithiophosphate Sodium Salt	A: 0.1M TEAA/ACN 99:1 B: 0.1M TEAA/ACN/MeOH 50:25:25	DVB/RP
Pg 93	Dilaurylthiodipropionate, Dimyristylthiodipropionate & Distearylthiodipropionate-polymer Antioxidants	45/45/10 ACN/IPA/IO	DVB RP 103Å
Pg 70	Dimethylaminoethyl Phenol Standards	80/20 ACN/0.2M NaOH	DVB/RP
Pg 134	Divinylbenzene	85/15 CHCl <sub>3</sub> /MeOH	DVB Sulfonated
Pg 30	Erythromycin (New & Aged)	45/23/32 ACN/0.05M KH <sub>2</sub> PO <sub>4</sub> /MeOH	DVB C18
Pg 7	Fat Soluble Vitamins	ACN w/0.1% TFA	DVB C18
Pg 66	Fatty Acids (Even Carbon)	95/5 MeOH/THF	DVB/RP
Pg 64	Fatty Acids (Odd Carbon)	95/5 MeOH/THF	DVB/RP
Pg 67	Short Chain Fatty Acids	50/50/1 ACN/H <sub>2</sub> O/H <sub>3</sub> PO <sub>4</sub>	DVB/RP
Pg 28	Short Chain Fatty Acids	50/50/1 ACN/H <sub>2</sub> O/H <sub>3</sub> PO <sub>4</sub>	DVB C18
Pg 76	Fumonisin B1	50/25/25 H <sub>2</sub> O/MeOH/ACN	DVB/RP
Pg 33	Fumonisin B1	50/25/25 H <sub>2</sub> O/MeOH/ACN	DVB C18
Pg 34	Fumonisin B1	50/25/25 H <sub>2</sub> O/MeOH/ACN	DVB C18
Pg 42	Fumonisin B1	50/25/25 H <sub>2</sub> O/MeOH/ACN	DVB C18
Pg 35	Fumonisin B1	50/25/25 H <sub>2</sub> O/MeOH/ACN	DVB C18
Pg 136	Glucose & Acetic Acid	0.01M Phosphoric Acid, pH adj. to 7.5	DVB Sulfonated
Pg 72	Glycerol Monoesters	78/22 to 100/0 ACN/H <sub>2</sub> O, 30 min.	DVB/RP
Pg 73	Glycerol Monoesters	78/22 to 100/0 ACN/H <sub>2</sub> O, 30 min.	DVB/RP
Pg 56	Hindered Phenolic Antioxidants	ACN/H <sub>2</sub> O 78/22 to 100/0 30 min linear	DVB/RP
Pg 98	Hydroquinone	10/10/80 MeOH/ACN/H <sub>2</sub> O w/0.1% TFA	DVB 103Å
Pg 94	Irganox 1010, Naugard 445 & Distearylthiodipropionate-Polymer Antioxidants	45/45/10 ACN/IPA/IO	DVB RP 103Å
Pg 121	Jordi Peptide Protein	90/10→60/40 A/B 30 min. linear A: 99/1 H <sub>2</sub> O/Acetic Acid B: 99/1 ACN/Acetic Acid	Jordi Peptide Protein
Pg 108	Lipids and a Triglyceride	50/43/7 CHCl <sub>3</sub> /MeOH/0.15% TFA in H <sub>2</sub> O	DVB Glucose NP 500Å
Pg 101	Mineral Oil and Motor Oils	75/25 Ethanol/2-propanol	DVB RP 103Å
Pg 22	2' and 3' Monophosphate Nucleotides	99/1 0.01M Sodium Acetate/ACN	DVB C18
Pg 23	2' and 3' Monophosphate Nucleotides	99/1 0.01M Sodium Acetate/ACN	DVB C18
Pg 14	Nucleosides and Bases	99/1 0.01M Sodium Acetate/ACN	DVB C18
Pg 15	Nucleosides and Bases	99/1 0.01M Sodium Acetate/ACN	DVB C18
Pg 16	Nucleosides and Bases	99/1 0.01M Sodium Acetate/ACN	DVB C18
Pg 63	Nucleosides and Bases	97/1/2 0.01M NaOAc/ACN/MeOH	DVB/RP
Pg 62	Nucleotides	97/1/2 0.01M NaOAc/ACN/MeOH	DVB/RP
Pg 102	Nylon 6 Monomer	10/90→100/0 CH <sub>3</sub> CN/H <sub>2</sub> O	DVB RP 103Å
Pg 50	Organic Acids	90/10 MeOH/THF	DVB/RP
Pg 51	Organic Acids	34.7/15.3/50 ACN/THF/H <sub>2</sub> O	DVB/RP
Pg 79	Parabens by GPC & LC	Acetonitrile, Methanol, THF or 50/50 Acetonitrile/Methanol	DVB 500Å
Pg 92	PEG 425	60/40 MeOH/H <sub>2</sub> O	DVB RP 103Å
Pg 31	Penicillin (New & Aged)	50/50 ACN/0.05M KH <sub>2</sub> PO <sub>4</sub>	DVB C18

Pg 59 Phenols	A: 45/45/10 B: 8/73/19 H <sub>2</sub> O/ACN/MeOH w/0.1% TFA	DVB/RP
Pg 24 Phenols	A: 44/56 B: 10/90 H <sub>2</sub> O/ACN w/0.1% TFA	DVB C18
Pg 129 Peptide Protein Standard Separation	80/20→40/60 A/B 30 min linear A: 0.15% TFA in H <sub>2</sub> O B: 0.15% TFA in ACN	Peptide/Protein Hydroxylated
Pg 107 Phospholipid Mixture	48/42/10 CHCl <sub>3</sub> /MeOH/H <sub>2</sub> O w/0.2% NH <sub>4</sub> OH	DVB Glucose 500Å
Pg 135 Phosphoric Compounds	0.02M Phosphoric Acid w/3% ACN	DVB Sulfonated
Pg 116 Polyhydroxy Alkaloids	A: ACN: 0.2% (w/v) B: Ammonium Acetate in H <sub>2</sub> O	Polyamine 500Å
Pg 99 Polymer Sample Separations	100/0→0/100 ACN/CHCl <sub>3</sub> linear over 30 minutes	DVB 103Å
Pg 105 Proteins	A: 0.1% TFA in H <sub>2</sub> O B: 0.1% TFA in ACN	DVB Solid Bead
Pg 104 Proteins & Enzymes	A: 0.1% TFA in H <sub>2</sub> O B: 0.1% TFA in ACN	DVB Solid Bead
Pg 131 Protein Standards Mix	80/20→40/60 A/B 30 min. Linear A: 0.15% TFA in H <sub>2</sub> O B: 0.15% TFA in ACN	Peptide Protein Hydroxylated
Pg 122 Protein Standard Separation	90/10→60/40 A/B 30 min. Linear A: 98/2 H <sub>2</sub> O/Acetic Acid B: 98/2 ACN/Acetic Acid	Peptide Protein
Pg 123 Protein Standard Separation	80/20→40/60 A/B 10 min. Linear A: 0.15% TFA in H <sub>2</sub> O B: 0.15% TFA in ACN	Peptide Protein
Pg 124 Protein Standard Separation	80/20→40/60 A/B 10 min. Linear A: 0.15% TFA in H <sub>2</sub> O B: 0.15% TFA in ACN	Peptide Protein
Pg 125 Protein Standard Separation	80/20→40/60 A/B 10 min. Linear A: 0.15% TFA in H <sub>2</sub> O B: 0.15% TFA in ACN	Peptide Protein
Pg 126 Protein Standard Separation	80/20→40/60 A/B 10 min. Linear A: 0.15% TFA in H <sub>2</sub> O B: 0.15% TFA in ACN	Peptide Protein
Pg 127 Protein Standard Separation	80/20→40/60 A/B 10 min. Linear A: 0.15% TFA in H <sub>2</sub> O B: 0.15% TFA in ACN	Peptide Protein
Pg 128 Protein Standard Separation	80/20→40/60 A/B 10 min. Linear A: 0.15% TFA in H <sub>2</sub> O B: 0.15% TFA in ACN	Peptide Protein
Pg 120 Protein Standard Separation	80/20→40/60 A/B 30 min. Linear A: 0.15% TFA in H <sub>2</sub> O B: 0.15% TFA in ACN	Peptide Protein
Pg 132 Protein Standard Separation	80/20→40/60 A/B 30 min. Linear A: 0.15% TFA in H <sub>2</sub> O B: 0.15% TFA in ACN	Peptide Protein
Pg 80 Quercitin & Rutin	75/25 Acetic Acid/MeOH	DVB 500Å
Pg 81 Quercitin & Rutin	50/50 Pyridine/Glacial Acetic Acid	DVB 500Å
Pg 95 RP Mix	MeOH/ACN/THF/H <sub>2</sub> O 10/20/30/40	DVB RP 103Å
Pg 41 RP Mix	MeOH/ACN/THF/H <sub>2</sub> O 10/20/30/40	DVB C18
Pg 88 RP Mix	MeOH/ACN/THF/H <sub>2</sub> O 10/20/30/40	DVB RP 500Å Bullet
Pg 38 RP Mix	MeOH/ACN/THF/H <sub>2</sub> O 10/20/30/40	DVB C18
Pg 77 RP Mix w/New Mobile Phase	50/25/25 ACN/MeOH/H <sub>2</sub> O	DVB/RP
Pg 78 RP Mix w/New Mobile Phase	30/30/40 THF/ACN/H <sub>2</sub> O	DVB/RP
Pg 40 RP Mix w/New Mobile Phase	30/30/40 THF/ACN/H <sub>2</sub> O	DVB C18

Pg 43	RP Mix w/New Mobile Phase	30/30/40 THF/ACN/H2O	DVB C18
Pg 36	RP Mix w/New Mobile Phase	60/40 THF/H2O	DVB C18
Pg 39	RP Mix w/New Mobile Phase	50/25/25 ACN/MeOH/H2O	DVB C18
Pg 96	Sea Nine 211 Rohm & Haas	100/0 CH3CN/H2O Antifouling Agent	DVB RP 103Å
Pg 112	Separation of Organic Acids	0.01M of Phosphoric Acid, pH 3 w/NaOH	DVB Organic Acid 500Å
Pg 89	Separation of Polyethylene Glycol & Polypropylene Glycol	A: 78/22 CH3Ca/H2O B: 2-propanol 40/60 A/B→B linear over 10 minutes	DVB 500Å
Pg 90	Separation of Polyethylene Glycol & Polypropylene Glycol	2-propanol	DVB 500Å
Pg 74	Slip Agents	78/22 to 100/0 ACN/H2O 30 min	DVB/RP
Pg 75	Slip Agents	78/22 to 100/0 ACN/H2O 30 min	DVB/RP
Pg 27	Steroids	65/35 ACN/H2O	DVB C18
Pg 65	Steroids	85/15 ACN/H2O	DVB/RP
Pg 117	Sugar Standards	75/20/5 ACN/H2O/MeOH	DVB Polyamine
Pg 114	Sugar Standards	80/15/5 ACN/H2O/MeOH	DVB Polyamine
Pg 71	Tannic Acid	80/10/10 H2O/ACN/MeOH	DVB/RP
Pg 68	Tetracycline & related Compounds	30/70 ACN/H2O w/0.1% TFA	DVB/RP
Pg 29	Tetracycline & related Compounds	20/80 ACN/H2O w/0.1% TFA	DVB C18
Pg 69	2,4,6-Tris(Dimethylaminoethyl)	80/20 ACN/0.02M NaOH Phenol Standard	DVB/RP
Pg 97	Urethane Prepolymer	100/0→0/100 ACN/CHCl3 linear over 30 min	DVB/103Å
Pg 55	Vanillin Compounds	80/20 ACN/H2O w/0.1% TFA	DVB/RP
Pg 13	Vanillin Compounds	40/60 H2O/ACN w/0.1% TFA	DVB C18
Pg 20	Vanillin Compounds	40/60 H2O/ACN w/0.1% TFA	DVB C18
Pg 12	Vasoconstrictors	24/75/1 ACN/0.2M NaOH/Butylamine	DVB C18
Pg 52	Vasoconstrictors	24/74/2 ACN/0.2M NaOH/Butylamine	DVB/RP
Pg 3	Water Soluble Vitamins	50/40/10 H2O/ACN/MeOH	DVB C18
Pg 4	Water Soluble Vitamins	50/40/10 H2O/ACN/MeOH	DVB C18
Pg 5	Water Soluble Vitamins	70/30 H2O/ACN w/0.1% TFA	DVB C18
Pg 6	Water Soluble Vitamins	50/40/10 H2O/ACN/MeOH	DVB C18