HPLC Detectors

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The Evaporative Lightscattering Detector ZAM 4000 is an universal detector for HPLC. It is used to analyze components which do not have UV absorption, cannot be separated with an isocratic solvent and must use a gradient elution, which cannot be used with a refractive index detector. Only mobile phases without any buffers should be used.

In contrast to other evaporative light scattering detectors all parts of the the ZAM 4000 detectors which come in contact with the sample are made of PTFE or glass. Therefore this parts are inert to common solvents and easy to clean. One more advantage of glass parts is the fact, that contamination can be located easily.

The nebulization takes place in the glass nebulization chamber. From here the aerosol is guided into the evaporation chamber which is made of glass, too. Here we add a **fixed, heated additiv gas flow of 1,9 l/min.** to the aerosol.

The solvent evaporated and only small sample particles reach the detection chamber. Here this particles pass a beam of bright light. By this the light is scattered in all directions. The scattered light is detected by a photomultiplier in an angel of 120°. The signal is amplified and exported as the detector signal.

All parts which get in contact with the sample can be unmounted and cleaned easily by user.

Specifications ZAM 4000

Detection Method :	Highly sensitive photomultiplier tube
Nebulizer Material Drying :	Glass, easy to clean, resistant to all common solvents
Chamber Material:	Glass, easy to clean, resistant to all common solvents
Nebulizer Chamber :	Glass, easy to clean, resistant to all common solvents
Material Light Source:	White LED
HPLC Mobile Phase Flow Rate:	0,5 – 3,5 ml/min.
Additiv Gas Controlling:	fix gas flow (1,9l/min.)
Gas consumption :	~ 3,7 l/min.
Recorder Output :	0 V - +1.2 V
Digital Interface:	RS232 serial port
Digital Input:	AutoZero, Start (Gas on / PMT power) (TTL)
Temperature Settings:	35°C - 85°C, ambient Temperature
Data Smoothing:	Raw, Fast, Slow, Medium
Power Source:	110 V AC / 220 V AC
Dimensions:	18cm x 49cm x50cm (W x H x D)
Weight:	24 kg