

UV/Vis - Detector S 4245 and S 4250

The UV/Vis - Detektors S 4245 und S 4250 are variable wavelength UV/Vis detectors for routine analysis and sophisticated research. The dual lamp design offers a wavelength range of 190 – 900 nm with a low baseline noise. The front - accessible flow cell can easily be exchanged, as can be the lamps which are accessible through a side panel in the instrument housing.



Integrated Wavelength Program

The UV/Vis - Detektoren S 4245 und S 4250 features a wavelength program to change the selected wavelength over time. With this feature the optimum wavelength can be selected for each analyzed substance according to its retention time.

Integrated Peak Detector

The integrated Peak Detector works as a basic fraction collector. The peak detection level can be freely programmed for peak start and peak end to enhance the collection purity. An integrated 24V output for switching a solenoid valve is used for the fraction collection, which is automatically operated with a selectable time delay.

Optional - Dual - Wavelength

The UV/Vis - Detektors S 4245 und S 4250 are available with an optional second wavelength. This feature enhances the Wavelength Program feature that you can measure 2 different wavelengths at the same time. A second D/A converter output comes with this option to keep the system flexible to be used with any data acquisition software available.

Optional - ONLINE - SCAN

Another option for the UV/Vis - Detektors S 4245 und S 4250 is the Online Scan. With the Online Scan whole spectrum information can be gathered at a certain time. This scan information is stored internally and can be accessed at any time. The Online Scan is a good alternative to a full UV PDA detector.




Figure: S 4245 / S4250Peak detection




Figure: Optical module S 4245 / S 4250

Technical Specifications*

| | |
|-----------------------------|--|
| Wetted Materials: | Stainless Steel / PEEK* |
| Baseline Noise: | $\pm 1 \times 10^{-5}$ AU (@240 nm, 1 sec. Rise time) |
| Baseline Drift: | 2×10^{-4} AU/h |
| Wavelength Range: | 190 – 900 nm |
| Wavelength Accuracy: | ± 2 nm |
| Linearity | > 2.0 AU |
| Light source: | Deuterium Lamp, Tungsten Lamp |
| Flow cell volume: | Micro 0,6 μ l, 3,0mm flow path Analytical 7,9 μ l, 10,0mm flow path Preparative 2,7 μ l, 1,5mm flow path |
| Wavelength Program: | Programmable, 10 Steps |
| Analog Output: | 1 x 1 V (optional: 2 x 1V) |
| Control Features: | Internal Peak Detector with +24 V solenoid switching output |
| Dimensions: | S 4245 310 x 165 x 478 mm (B x H x T) S 4250 396 x 165 x 478 mm (B x H x T) |
| Power Supply: | 100 - 250 V (47—63 Hz) |

*abhängig von der Gerätekonfiguration

Order Information UV/Vis - Detector S 4245 and S 4250

| Part-No.: | Description: |
|----------------|--|
| UV/Vis 4245 | UV/Vis Detector S 4245, 1 Channel |
| UV/Vis 4245-2 | UV/Vis Detector S 4245, 2 Channels |
| UV/Vis 4245-2S | UV/Vis Detector S 4245, 2 Channels with SCAN Option |
| UV/Vis 4250 | UV/Vis Detector S 4250, 1 Channel |
| UV/Vis 4250-2 | UV/Vis Detector S 4250, 2 Channel |
| UV/Vis 4250-2S | UV/Vis Detector S 4250, 2 Channels with SCAN Option |
| UV/VIS—1032009 | S 4245 / S 4250 Flow cell, analytical, stainless Steel |
| UV/VIS—1032010 | S 4245 / S 4250 Flow cell, analytical, PEEK |
| UV/VIS—1032011 | S 4245 / S 4250 Flow cell, micro, stainless Steel |
| UV/VIS—1032012 | S 4245 / S 4250 Flow cell, micro, PEEK |
| UV/VIS—1032013 | S 4245 / S 4250 Flow cell, semi-preparative, stainless Steel |
| UV/VIS—1032014 | S 4250 / S 4250 Flow cell, semi-preparative, stainless Steel |

Order Information Spare Parts UV/Vis - Detector S 4245 and S 4250

| Part-No.: | Description: , |
|-------------|---|
| 47- 4010001 | Tungsten Lamp for S 4245 / S 4250, preadjusted |
| 47- 4010002 | Deuterium Lamp for S 4245 / S 4250, preadjusted |