# HPLC Detectors

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The Evaporative Lightscattering Detector ZAM 3000 is an universal detector for HPLC. It is used to analyze components which do not have UV absorption, cannot be separated with an isocratic solvent and must use a gradient elution, which cannot be used with a refractive index detector. Only mobile phases without any buffers should be used.

In contrast to other evaporative light scattering detectors all parts of the the ZAM 3000 detectors which come in contact with the sample are made of PTFE or glass. Therefore this parts are inert to common solvents and easy to clean. One more advantage of glass parts is the fact, that contamination can be located easily.

The nebulization takes place in the glass nebulization chamber. From here the aerosol is guided into the evaporation chamber which is made of glass, too. **Here we add, with a mass flow controller heated additiv gas to the aerosol. The heated additiv gas is variable adjustable (0,5 - 3,5l/min.).**

The solvent evaporated and only small sample particles reach the detection chamber. Here this particles pass a beam of bright light. By this the light is scattered in all directions. The scattered light is detected by a photomultiplier in an angel of 120°. The signal is amplified and exported as the detector signal.

All parts which get in contact with the sample can be unmounted and cleaned easily by user.

# Specifications ZAM3000

<b>Detection Method :</b>	Highly sensitive photomultiplier tube
<b>Nebulizer Material Drying :</b>	Glass, easy to clean, resistant to all common solvents
<b>Chamber Material:</b>	Glass, easy to clean, resistant to all common solvents
<b>Nebulizer Chamber :</b>	Glass, easy to clean, resistant to all common solvents
<b>Material Light Source:</b>	White LED
<b>HPLC Mobile Phase Flow Rate:</b>	0,5 – 3,5 ml/min.
<b>Additiv Gas Controlling:</b>	by Mass flow controller (0,5 – 3,5 l/min.)
<b>Gas Consumption :</b>	~ 1,7 – 4,5 l/min.
<b>Recorder Output :</b>	0 V - +1.2 V
<b>Digital Interface:</b>	RS232 serial port
<b>Digital Input:</b>	AutoZero, Start (Gas on / PMT power) (TTL)
<b>Temperature Settings:</b>	35°C - 85°C, ambient Temperature
<b>Data Smoothing:</b>	Raw, Fast, Slow, Medium
<b>Power Source:</b>	110 V AC / 220 V AC
<b>Dimensions:</b>	18cm W / 49cm H / 50cm D
<b>Weight:</b>	24 kg