



SALE ! **Demo unit YL9120 UV/Vis - Detector** **SALE !**
including YL-Clarity data acquisition system (unit run time < 30 hours)



Wavelength Range :	190-900 nm
Data collection rate :	up to 50Hz
Light Source :	Deuterium lamp & tungsten lamp
Noise level :	< 0.5 X 10 ⁻⁵ AU , 254nm, dry cell
Drift :	<1 X 10 ⁻⁴ AU/hr
Bandwidth :	5.5 nm
Wavelength Accuracy :	1 nm
Wavelength Precision :	0.1nm
Linearity :	> 99.5% for 2.5 AU (acetone, 254nm)
Path Length :	10 mm (Analytical cell)
Cell Volume :	10 uL (Analytical cell)

Special Price for Demo Unit: 3999,00 €, net, ex work, full warranty

List price: 5.400,00€ you save more than 25%

Schambeck SFD GmbH, Drieschweg 13A, D - 53604 Bad Honnef, Phone: ++49-2224-9239-0,
Fax: ++49-2224-9239-20, e-mail: mailto@schambeck-sfd.com, Internet: <http://www.schambeck-sfd.com>
Geschäftsführer: Karl-Heinz Schambeck, Amtsgericht Siegburg: HRB 6811, VAT: DE 122267991
Bankverbindung: Stadtparkasse Bad Honnef, BLZ 380 512 90, Konto-Nr. 153 320